


## Anti-Phospho-RAD51 (Thr309) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Applications

Verified Activity: 1. Western blot analysis of HELA-UV 293T-UV using p-Rad51 (T309) antibody.



Application:	ELISA,IF,IHC,WB
Recommended	IF: 1:50-200; WB: 1:500-2000; IHC: 1:100-300; ELISA: 1:5000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthesized phosphopeptide: human Rad51 around the phosphorylation site of T309
Antigen Species:	Human
Uniprot ID:	Q06609
Synonyms:	RAD51 (p-T309);RAD51 (p-Thr309);p-RAD51 (Thr309);p-RAD51 (T309)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481