


## Anti-Phospho-BRAF (Ser602) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 86 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Applications

Verified Activity: 1. Western Blot analysis of 3T3 cells using Phospho-Raf-B (S602) Polyclonal Antibody.



Application: ELISA,IHC,WB  
Recommended WB: 1:500-2000; IHC: 1:100-300; ELISA: 1:20000

## Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

## Antigen Details

Immunogen: A synthesized phosphopeptide: human Raf-B around the phosphorylation site of S602  
Antigen Species: Human  
Uniprot ID: P15056  
Synonyms: p-BRAF (S602);BRAF (p-Ser602);p-BRAF (Ser602);BRAF (p-S602)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481