

Anti-Flag Tag Antibody-HRP (2U342)

Product Details

Reactivity:	Species independent
Conjugation:	HRP
Clone:	2U342
Purification:	affinity purified by Protein A

Applications

Application:	ELISA,WB
Recommended	WB: 1:30000-100000; ELISA: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481