

## Anti-HA Tag Antibody-HRP (4M675)

### Product Details

Reactivity:	Species independent
Conjugation:	HRP
Clone:	4M675
Purification:	affinity purified by Protein A

### Applications

Application:	ELISA,WB
Recommended	WB: 1:30000-100000; ELISA: 1:5000-10000

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. Keep away from direct sunlight.
Shipping:	Shipping with blue ice.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481