


Anti-RFP Polyclonal Antibody

Product Details

Conjugation: Unconjugated
Purification: Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of RFP Recombinant protein, using TMAC-03590 diluted at 1) 1:5,000, 2) 1:10,000.



Application: IHC,IP,WB
Recommended WB: 1:2000-5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

Research Background

The red fluorescent protein cloned from *Discosoma coral* (DsRed or drFP583) holds great promise for biotechnology and cell biology as a spectrally distinct companion or substitute for the green fluorescent protein (GFP) from the *Aequorea jellyfish*.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481