

Anti-Phospho-Rrn3 (Ser649) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 70 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Application:	ELISA,IF,IHC-P,WB
Recommended	WB: 1:500-2000; IHC-P: 1:100-300; ELISA: 1:10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human TIF-IA around the phosphorylation site of Ser649. AA range:602-651
Antigen Species:	Human
Uniprot ID:	Q9NYV6
Synonyms:	Rrn3 (p-Ser649);p-Rrn3 (S649);Rrn3 (p-S649);p-Rrn3 (Ser649)

Research Background

function:Required for efficient transcription initiation by RNA polymerase I.,similarity:Belongs to the RRN3 family.,

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481