

Anti-SMCC Antibody (8U377)

Product Details

Reactivity:	Species independent
Conjugation:	Unconjugated
Clone:	8U377
Purification:	ProA affinity purified

Applications

Application:	ELISA
Recommended	ELISA: 1:50000-100000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months (Store at 4°C after thawing). Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein: SMCC
Synonyms:	22 kDa actin binding protein;SM22;Transgelin variant 2;SM22-alpha;Human 22kDa smooth muscle protein;TAGLN;Transgelin;22 kDa actin-binding protein;WS3 10;TAGL_HUMAN;SMCC;SM22 alpha;Protein WS3-10;TAGLN1;Smooth muscle protein 22-alpha

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481