


## Anti-SRT Tag Antibody (6Q868)

### Product Details

Conjugation: Unconjugated  
Clone: 6Q868  
Purification: Affinity purification using immunogen.

### Applications

Verified Activity: 1. Western blot analysis of SRT-Recombinant protein, using TMAC-03858 diluted at 1) 1:5,000 2) 1:10,000.



Application: WB  
Recommended WB: 1:5000-10000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Research Background

The SRT-tag is a peptide epitope comprised of 10 amino acids (TFIGAIATDT) derived from the crystalline surface layer protein of Rickettsia typhi1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481