



## Anti-Survivin Antibody (9D807)

### Product Details

Reactivity:	Human,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 16 kDa.
Clone:	9D807
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

### Applications

- Verified Activity:
1. Western blot analysis of 1) HeLa, 2) 293, 3) PC12 using Survivin Monoclonal Antibody.
  2. Immunohistochemical analysis of paraffin-embedded human colon carcinoma using Survivin Monoclonal Antibody.



Application:	IHC,WB
Recommended	WB: 1:1000-2000; IHC: 1:200-500

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein

Uniprot ID: O15392

Synonyms: IAP4;baculoviral IAP repeat containing 5;API4;BIRC5;EPR-1

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481