


## Anti-Phospho-Syk (Tyr525) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat,Monkey
Conjugation:	Unconjugated
Molecular Weight:	Actual: 72 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Applications

Verified Activity: 1. Western Blot analysis of COS7 cells using Phospho-Syk (Y525) Polyclonal Antibody.



Application:	ELISA,IHC,WB
Recommended	WB: 1:500-2000; IHC: 1:100-300; ELISA: 1:20000

## Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

## Antigen Details

Immunogen:	A synthesized phosphopeptide: human Syk around the phosphorylation site of Y525
Antigen Species:	Human
Uniprot ID:	P43405
Synonyms:	Tyrosine protein kinase SYK;40 kDa tyrosine-protein kinase SYK;Syk;Tyrosine-protein kinase SYK;p72syk;KSYK;p-Syk (Tyr525);Syk (p-Tyr525);Syk (p-Y525);p-Syk (Y525);Tyrosine-protein kinase SYK(Spleen tyrosine kinase) [Contains: 72kDa tyrosine-protein kinase SYK;Spleen tyrosine kinase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481