


## Anti-T7 Tag Antibody (5K790)

### Product Details

Ig Type:	IgG1
Conjugation:	Unconjugated
Clone:	5K790
Purification:	Affinity-chromatography

### Applications

Verified Activity: 1. Western blot analysis of GST-T7 fusion protein, using T7-tag mouse monoclonal antibody.



Application:	ICC/IF,WB
Recommended	WB: 1:3000-10000; ICC/IF: 1:200

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen:	A synthesized peptide: MASMTGGQQMG
Synonyms:	T 7 Tag

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481