

Anti-TAP Tag Antibody (2A205)

Product Details

Conjugation:	Unconjugated
Clone:	2A205
Purification:	Affinity purification using immunogen.

Applications

Verified Activity: 1. Western blot analysis of TAP recombinant protein, using TMAC-03965 diluted at 1:5,000.



Application:	IP,WB
Recommended	WB: 1:5000-10000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Research Background

Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481