

Anti-TBP Antibody (1C626)



Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 35-45 kDa.
Clone:	1C626

Applications

Verified Activity:

1. Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using TBP/TATA Binding Protein mAb diluted at 1:200.
2. Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using TBP/TATA Binding Protein Mouse mAb diluted at 1:200.



Application:	IHC
Recommended	IHC: 1:100-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: Recombinant Protein
Uniprot ID: P20226
Synonyms: TATA-box-binding protein;TATA-box factor;TFIID;TF2D;TATA sequence-binding protein;TATA-binding factor;GTF2D1;Transcription initiation factor TFIID TBP subunit

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481