

Anti-TBP Antibody (20811)


Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 35-45 kDa.
Clone:	20811

Applications

Verified Activity:

1. Western blot analysis of 1) HeLa, 2) Mouse Brain Tissue, 3) Rat Brain Tissue with TBP/TATA Binding Protein Mouse mAb diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:1000-2000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	P20226
Synonyms:	Transcription initiation factor TFIID TBP subunit;TATA sequence-binding protein;GTF2D1;TATA-box factor;TFIID;TATA-binding factor;TF2D;TATA-box-binding protein

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481