


## Anti-Phospho-TIRAP (Tyr86) Polyclonal Antibody

## Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 23 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

## Applications

Verified Activity: 1. Western blot analysis of K562 using p-TIRAP (Y86) antibody.



Application: ELISA,IHC,WB  
Recommended WB: 1:500-2000; IHC: 1:100-300;ELISA: 1:5000

## Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

## Antigen Details

Immunogen: A synthesized phosphopeptide: human TIRAP around the phosphorylation site of Y86  
Antigen Species: Human  
Uniprot ID: P58753  
Synonyms: TIRAP (p-Y86);TIRAP (p-Tyr86);p-TIRAP (Y86);p-TIRAP (Tyr86)

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481