

## Anti-TNF alpha Antibody (4B973)

### Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 17, 25 kDa.
Clone:	4B973

### Applications

Verified Activity:


1. Western blot analysis of Recombinant Human TNF a Protein with TNF  $\alpha$  Mouse mAb at 1:2,000.
2. The picture was kindly provided by our customer.
3. Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using TNF  $\alpha$  Mouse mAb diluted at 1:50.



mouse BV2



TNF  $\alpha$



Application:

IF,IHC,WB

Recommended

WB: 1:1000-3000; IHC: 1:50-200; IF: 1:50-200

**Properties**

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.

Shipping: Shipping with blue ice.

---

**Antigen Details**

Immunogen: A synthesized peptide

Uniprot ID: P01375

Synonyms: APC1 protein

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481