

Anti-TNF alpha Antibody (6X125)


Product Details

Ig Type:	IgG1
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 17, 25 kDa.
Clone:	6X125

Applications

Verified Activity:

1. Western blot analysis of Recombinant Human TNF α Protein with TNF α Mouse mAb diluted at 1:2,000.
2. Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using TNF α Mouse mAb diluted at 1:50.



Application:	IHC,WB
Recommended	WB: 1:1000-3000; IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide

Uniprot ID: P01375

Synonyms: APC1 protein

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481