

Anti-Transferrin Antibody (2K968)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 77 kDa.
Clone:	2K968

Applications

Verified Activity: Western blot analysis of Human Serum with Transferrin Mouse mAb diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:2000-5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	P02787
Synonyms:	transferrin;TFQTL1;PRO1557;HEL-S-71p;PRO2086

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481