

Anti-Transferrin Antibody (7Q166)

Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 77 kDa.
Clone:	7Q166

Applications

Verified Activity:	Immunohistochemical analysis of paraffin-embedded rat-liver using antibody diluted at 1:50.
Application:	IHC
Recommended	IHC: 1:100-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	P02787
Synonyms:	PRO2086;transferrin;HEL-S-71p;TFQTL1;PRO1557

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481