


## Anti-Transferrin Antibody (9Z409)

### Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 77 kDa.
Clone:	9Z409

### Applications

Verified Activity: 1. Western blot analysis of Human Serum with Transferrin Mouse mAb diluted at 1:2,000.



Application:	WB
Recommended	WB: 1:2000-5000

### Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.  
Shipping: Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	P02787
Synonyms:	PRO2086;transferrin;HEL-S-71p;TFQTL1;PRO1557

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481