

## Anti-Transferrin Antibody (7V26)

### Product Details

Ig Type:	IgG1
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 77 kDa.
Clone:	7V26

### Applications

Verified Activity:	1. Immunohistochemical analysis of paraffin-embedded human-liver using antibody diluted at 1:50. 2. IHC: 1:100-1:200.
Application:	IHC
Recommended	IHC: 1:100-200

### Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

### Antigen Details

Immunogen:	Recombinant Protein
Uniprot ID:	P02787
Synonyms:	HEL-S-71p;PRO1557;transferrin;TFQTL1;PRO2086

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481