


Anti-Phospho-NTRK2 (Tyr516) Polyclonal Antibody

Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 92 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity: 1. Western Blot analysis of RAT-MUSCLE cells using Phospho-Trk B (Y516) Polyclonal Antibody.



Application: ELISA,IHC,WB
Recommended WB: 1:500-2000; IHC: 1:100-300; ELISA: 1:5000

Properties

Stability & Storage: Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping: Shipping with blue ice.

Antigen Details

Immunogen: A synthesized phosphopeptide: human Trk B around the phosphorylation site of Y516
Antigen Species: human
Uniprot ID: Q16620
Synonyms: NTRK2 (p-Tyr516);NTRK2 (p-Y516);p-NTRK2 (Tyr516);p-NTRK2 (Y516)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481