

Anti-Tyrosine Hydroxylase Antibody (7C683)

Product Details

Ig Type:	IgG1,Kappa
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 58 kDa.
Clone:	7C683
Purification:	Protein G

Applications

Application:	ELISA,IF,WB
Recommended	WB: 1:500-2000; IF: 1:100-500; ELISA 1:1000-5000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human Tyrosine Hydrolase
Antigen Species:	human
Uniprot ID:	P07101
Synonyms:	tyrosine hydroxylase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481