

Anti-Acetyl-UBB (Lys27) Polyclonal Antibody



Product Details

Ig Type:	IgG
Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 80 50 kDa.
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.

Applications

Verified Activity:

1. Western Blot analysis of HEPG2 cells using Acetyl-Ub (K27) Polyclonal Antibody.
2. Western Blot analysis of HeLa cells using Acetyl-Ub (K27) Polyclonal Antibody.



Application:	ELISA,WB
Recommended	WB: 1:500-2000; ELISA: 1:20000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized acetylpeptide: human Ub around the acetylation site of K27
Antigen Species: Human
Synonyms: Acetyl-UBB (K27);Ac-UBBK27;Ac-UBB (K27);UBBK27ac;UBBK27-acetylated;Ac-UBB (Lys27)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481