

Anti-V5 tag Antibody (9D730)


Product Details

Ig Type:	IgG
Conjugation:	Unconjugated
Clone:	9D730
Purification:	ProA affinity purified

Applications

Verified Activity:

1. Western blot analysis of V5 tag on V5-tagged recombinant protein lysates using anti-V5 tag antibody at 1/1,000 dilution.



Application:	WB
Recommended	WB: 1:1000-5000; ICC: 1:100-500; FCM: 1:50-100

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	Recombinant Protein
Synonyms:	V 5 tag

Research Background

The V5 epitope tag is derived from a small epitope (Pk) present on the P and V proteins of the paramyxovirus of simian virus 5 (SV5). The V5 antibody detects recombinant proteins containing the 14 amino acid GKPIPPLLGLDST sequence.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481