


Anti-VEGFA Antibody (9D606)

Product Details

| | |
|-------------------|--|
| Ig Type: | Mouse IgG1, Kappa |
| Reactivity: | Human |
| Conjugation: | Unconjugated |
| Molecular Weight: | Theoretical: 21(monomer), 42(dimer) kDa. |
| Clone: | 9D606 |
| Purification: | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |

Applications

Verified Activity: 1. Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-VEGF-A antibody.



| | |
|--------------|-------------------------------------|
| Application: | ELISA,WB |
| Recommended | WB: 1:500-2000; ELISA: 1:5000-20000 |

Properties

| | |
|----------------------|---|
| Stability & Storage: | Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles. |
| Shipping: | Shipping with blue ice. |

Antigen Details

| | |
|------------------|-------------------------------------|
| Immunogen: | A synthesized peptide: human VEGF-A |
| Antigen Species: | Human |
| Gene ID: | 7422 |
| Uniprot ID: | P15692 |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481