


Anti-VEGFA Antibody (9D606)

Product Details

Ig Type:	Mouse IgG1, Kappa
Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 21(monomer), 42(dimer) kDa.
Clone:	9D606
Purification:	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.

Applications

Verified Activity: 1. Whole cell lysates were separated by 10% SDS-PAGE, and the membrane was blotted with anti-VEGF-A antibody.



Application:	ELISA,WB
Recommended	WB: 1:500-2000; ELISA: 1:5000-20000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized peptide: human VEGF-A
Antigen Species:	human
Gene ID:	7422
Uniprot ID:	P15692

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481