

Anti-XRCC4 Antibody (2Z65)



Product Details

Reactivity:	Human
Conjugation:	Unconjugated
Molecular Weight:	Actual: 38 kDa.
Clone:	2Z65
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

Verified Activity:

- 1) Input: HeLa Cell Lysate 2) IP product: IP dilute 1:200.
2. Western blot analysis of 1) HeLa, 2) 293T, diluted at 1:3000.



Application:	IP,WB
Recommended	WB: 1:2000; IP: 1:200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Uniprot ID: Q13426

Synonyms: X-ray repair cross-complementing protein 4;XRCC 4;DNA repair protein XRCC4

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481