



Anti-ACTA2 Antibody (2K214)

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Actual: 42 kDa.
Clone:	2K214
Purification:	The antibody was affinity-purified from antiserum by affinity-chromatography

Applications

- Verified Activity:
1. Western blot analysis of 1) Hela, 2) 3T3, 3) Rat Brain using α -SMA Monoclonal Antibody.
 2. Immunohistochemical analysis of paraffin-embedded Mouse Cecal Tissue using α -SMA Monoclonal Antibody.



Application:	IF,IHC,WB
Recommended	WB: 1:5000-50000; IHC: 1:1000-2000; IF: 1:200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide
Uniprot ID: P62736
Synonyms: ACTA2;Cell growth-inhibiting gene 46 protein;aortic smooth muscle;Actin, aortic smooth muscle;ACTSA;ACTVS;Alpha-actin-2

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481