

Anti-Beta-Actin Antibody (5H32)



Product Details

Ig Type:	IgG1
Reactivity:	Human, Mouse, Rat, Monkey, Dog, Chicken, Hamouseter, Rabbit, Insect
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 43 kDa.
Clone:	5H32

Applications

Verified Activity:

1. Western blot analysis of Hela (1), Rat brain (2), Mouse brain (3), Chicken lung (4), and Rabbit testis (5) diluted at 1:10000.
2. IHC Staining of Human ovary tissue, diluted at 1:200.



Application:	IHC, WB
Recommended	WB: 1:500-2000; IHC: 1:50-200

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen: A synthesized peptide
Uniprot ID: P60709
Synonyms: β -actin; β -Actin Protein;Actin Cytoplasmic 1;Beta-Actin;ACTB

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481