

Anti-Phospho-CTNNB1 (Ser33/37/Thr41) Polyclonal Antibody

Product Details

Reactivity:	Human,Mouse,Rat
Conjugation:	Unconjugated
Molecular Weight:	Theoretical: 92 kDa.
Purification:	Affinity purification using immunogen.

Applications

Application:	WB
Recommended	WB: 1:1000-3000

Properties

Stability & Storage:	Store at -20°C or -80°C for 12 months. Avoid repeated freeze-thaw cycles.
Shipping:	Shipping with blue ice.

Antigen Details

Immunogen:	A synthesized phosphopeptide: human β -Catenin around the phosphorylation site of Ser33/37/Thr41
Antigen Species:	human
Uniprot ID:	P35222
Synonyms:	CTNNB1 (p-S33/37/T41);CTNNB1 (p-Ser33/37/Thr41);p-CTNNB1 (Ser33/37/Thr41);p-CTNNB1 (S33/37/T41)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481