

(2E)-1-(2-hydroxyphenyl)-3-phenylprop-2-en-1-one

Chemical Properties


CAS No. :

Formula: C₁₅H₁₂O₂

Molecular Weight: 224.26

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	(2E)-1-(2-hydroxyphenyl)-3-phenylprop-2-en-1-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.4591 mL	22.2955 mL	44.5911 mL
5 mM	0.8918 mL	4.4591 mL	8.9182 mL
10 mM	0.4459 mL	2.2296 mL	4.4591 mL
50 mM	0.0892 mL	0.4459 mL	0.8918 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481