


3'-dC(Bz)-2'-phosphoramidite

Chemical Properties

CAS No. :	157327-96-3
Formula:	C ₄₆ H ₅₂ N ₅ O ₈ P
Molecular Weight:	833.91
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3'-dC(Bz)-2'-phosphoramidite is a Nucleoside Phosphoramidite.
Targets(IC50)	Nucleoside Antimetabolite/Analog,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1992 mL	5.9959 mL	11.9917 mL
5 mM	0.2398 mL	1.1992 mL	2.3983 mL
10 mM	0.1199 mL	0.5996 mL	1.1992 mL
50 mM	0.024 mL	0.1199 mL	0.2398 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481