

HNRNPL Protein, Human, Recombinant (aa 1-589, His & SUMO)

General Information

Synonyms: HNRNPL;HNRPL;Heterogeneous nuclear ribonucleoprotein L;hnRNP L

Protein Construction: 1-589 aa

Species: Human

Expression Host: E. coli

Accession: P14866

Molecular Weight: 77.1 kDa (predicted)

AA Sequence: MSRRLLPRAEKRRRRLERQRRRSGAMVKMAAAGGGGGGGRYGGGSEGGRAPKRLKTDNAGDQ
HGGGGGGGGGAGAAGGGGGGGENYDDPHKTPASPVVHIRGLIDGVVEADLVEALQEFGPISYVVVMPKKRQ
ALVEFEDVLGACNAVNYAADNQIYIAGHPAFVNYSTSQKISRPGSDSDSRVNSVLLFTILNPIYSITTDVLYTIC
NPCGPVQRIVIFRKNQVQAMVEFDSVQSAQRAKASLNGADIYSGCCTLKIEYAKPTRLNVFKNDQDTWDYTN
PNLSGQGDGPGSNPNKRQRQPPLGDHPAEYGGPHGGYHSHYHDEGYGPPPPHYEGRRMGPPVGGHRRGP
SRYGPQYGHPPPPPPPEYGPHADSPVLMVYGLDQSKMNCDFVNFVFLYGNVEKVKFMKSKPGAAMVEM
ADGYAVDRAITHLNNFMFGQKLNVCVSKQPAIMPGQSYGLEDGSCSYKDFSESRRNRFSTPEQAANKRIQ
HPSNVLHFFNAPLEVTEENFFEICDELGVKRPSSVKVFSGKSERSSSGLLEWESKSDALETGLFLNHYQMKNP
NGPYPYTLKLCFSTAQHAS

QC Testing

Biological Activity: Activity has not been tested. It is theoretically active, but we cannot guarantee it. If you require protein activity, we recommend choosing the eukaryotic expression version first.

Purity: > 90% as determined by SDS-PAGE.

Endotoxin: < 1.0 EU/μg of the protein as determined by the LAL method.

Formulation: Lyophilized from a 0.2 μm sterile filtered PBS, 6% Trehalose, pH 7.4.

Preparation and Storage

Reconstitution:

Reconstitute the lyophilized protein in distilled water. The product concentration should not be less than 100 μg/ml. Before opening, centrifuge the tube to collect powder at the bottom. After adding the reconstitution buffer, avoid vortexing or pipetting for mixing.

Stability & Storage:

It is recommended to store recombinant proteins at -20°C to -80°C for future use. Lyophilized powders can be stably stored for over 12 months, while liquid products can be stored for 6-12 months at -80°C. For reconstituted protein solutions, the solution can be stored at -20°C to -80°C for at least 3 months. Please avoid multiple freeze-thaw cycles and store products in aliquots.

Actual storage temperature shall be subject to the COA.

Shipping:

In general, lyophilized powders are shipped with blue ice, while solutions are shipped with dry ice.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481