


Senegose J

Chemical Properties

CAS No. :	156250-46-3
Formula:	C60H74O34
Molecular Weight:	1339.21
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Senegose J is an oligosaccharide component extracted from Polygala senega var. latifolia.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7467 mL	3.7335 mL	7.4671 mL
5 mM	0.1493 mL	0.7467 mL	1.4934 mL
10 mM	0.0747 mL	0.3734 mL	0.7467 mL
50 mM	0.0149 mL	0.0747 mL	0.1493 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481