


Bekanamycin sulfate

Chemical Properties

CAS No. : 29701-07-3
 Formula: C₁₈H₃₉N₅O₁₄S
 Molecular Weight: 581.59
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bekanamycin sulfate is an agent of bioactive chemicals.
Targets(IC50)	Others,Antibacterial,Antibiotic

Solubility Information

Solubility	DMSO: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7194 mL	8.5971 mL	17.1942 mL
5 mM	0.3439 mL	1.7194 mL	3.4388 mL
10 mM	0.1719 mL	0.8597 mL	1.7194 mL
50 mM	0.0344 mL	0.1719 mL	0.3439 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481