

Fmoc-PEG3-Ala-Ala-Asn(Trt)-PAB

Chemical Properties

CAS No. : 2055042-69-6


Formula: C60H66N6O11

Molecular Weight: 1047.2

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Fmoc-PEG3-Ala-Ala-Asn(Trt)-PAB is a cleavable three-unit polyethylene glycol (PEG) linker used in the synthesis of antibody-drug conjugates (ADCs)[1].
Targets(IC50)	ADC Linker,PROTAC Linker
In vitro	ADCs consist of an antibody connected to an ADC cytotoxin via an ADC linker[1].

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.9549 mL	4.7746 mL	9.5493 mL
5 mM	0.191 mL	0.9549 mL	1.9099 mL
10 mM	0.0955 mL	0.4775 mL	0.9549 mL
50 mM	0.0191 mL	0.0955 mL	0.191 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Beck A, et al. Strategies and challenges for the next generation of antibody-drug conjugates. Nat Rev Drug Discov. 2017 May;16(5):315-337.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481