

Compound T130886

Chemical Properties


CAS No. :

Formula: C₂₆H₂₈O₁₅

Molecular Weight: 580.5

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description

4-(3,4-dihydroxyphenyl)-7-hydroxy-5-[(3,4,5-trihydroxy-6-[(3,4,5-trihydroxyoxan-2-yl)oxy]methyl]oxan-2-yl]oxy]-2H-chromen-2-one is a high-purity biochemical reagent suitable for biochemical experiments and drug synthesis research.

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7227 mL	8.6133 mL	17.2265 mL
5 mM	0.3445 mL	1.7227 mL	3.4453 mL
10 mM	0.1723 mL	0.8613 mL	1.7227 mL
50 mM	0.0345 mL	0.1723 mL	0.3445 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481