

## Epimedin A1

### Chemical Properties

CAS No. : 140147-77-9


Formula: C<sub>39</sub>H<sub>50</sub>O<sub>20</sub>

Molecular Weight: 838.80

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



### Biological Description

Description	Epimedin A1 (Hexandraside F)(Hexandraside F) is a flavonoid extracted from Herba Epimedii which is one of commonly used Chinese medicines.
Targets(IC50)	Others

### Solubility Information

Solubility	DMSO: 140.00 mg/mL (166.91 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2.00 mg/mL (2.38 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

### Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.1922 mL	5.9609 mL	11.9218 mL
5 mM	0.2384 mL	1.1922 mL	2.3844 mL
10 mM	0.1192 mL	0.5961 mL	1.1922 mL
50 mM	0.0238 mL	0.1192 mL	0.2384 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

### Reference

- Chen XJ, et al. A rapid method for simultaneous determination of 15 flavonoids in Epimedium using pressurized liquid extraction and ultra-performance liquid chromatography. *J Pharm Biomed Anal.* 2008 Jan 22;46(2):226-35.
- Gao Y, Xu G, Ma L, et al. Icariside I specifically facilitates ATP or nigericin-induced NLRP3 inflammasome activation and causes idiosyncratic hepatotoxicity. *Cell Communication and Signaling.* 2021 Feb 11;19(1):13. doi: 10.1186/s12964-020-00647-1.
- Gao Y, Xu G, Ma L, et al. Icarisid I specifically facilitates ATP or nigericin-induced NLRP3 inflammasome activation and causes idiosyncratic hepatotoxicity. *Cell Communication and Signaling.* 2020
- Gao Y, Xu G, Ma L, et al. Icarisid I specifically facilitates ATP or nigericin-induced NLRP3 inflammasome activation and causes idiosyncratic hepatotoxicity[J]. *Cell Communication and Signaling.* 2020
- Gao Y, Xu G, Ma L, et al. Icariside I specifically facilitates ATP or nigericin-induced NLRP3 inflammasome activation and causes idiosyncratic hepatotoxicity[J]. *Cell Communication and Signaling.* 2021, 19(1): 1-14.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481