


Metralindole

Chemical Properties

CAS No. :	54188-38-4
Formula:	C ₁₅ H ₁₇ N ₃ O
Molecular Weight:	255.32
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Metralindole is a reversible inhibitor of monoamine oxidase A (RIMA).
Targets(IC50)	Others,CDK

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.9167 mL	19.5833 mL	39.1665 mL
5 mM	0.7833 mL	3.9167 mL	7.8333 mL
10 mM	0.3917 mL	1.9583 mL	3.9167 mL
50 mM	0.0783 mL	0.3917 mL	0.7833 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481