


Benzenesulfonic acid, m-iodo-, sodium salt

Chemical Properties

CAS No. :	51119-76-7
Formula:	C ₆ H ₄ INaO ₃ S
Molecular Weight:	306.05
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Benzenesulfonic acid, m-iodo-, sodium salt is a bioactive chemical.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2674 mL	16.3372 mL	32.6744 mL
5 mM	0.6535 mL	3.2674 mL	6.5349 mL
10 mM	0.3267 mL	1.6337 mL	3.2674 mL
50 mM	0.0653 mL	0.3267 mL	0.6535 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481