

## Isamoltan

## Chemical Properties


CAS No. : 116861-00-8

Formula: C<sub>16</sub>H<sub>22</sub>N<sub>2</sub>O<sub>2</sub>

Molecular Weight: 274.36

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Isamoltan is a 5-HT(1B) receptor antagonist.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.6448 mL	18.2242 mL	36.4485 mL
5 mM	0.729 mL	3.6448 mL	7.2897 mL
10 mM	0.3645 mL	1.8224 mL	3.6448 mL
50 mM	0.0729 mL	0.3645 mL	0.729 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481