


Isobutyl Salicylate

Chemical Properties

CAS No. :	87-19-4
Formula:	C ₁₁ H ₁₄ O ₃
Molecular Weight:	194.23
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Isobutyl Salicylate is a phenolic benzoate ester functionally related to salicylic acid. Isobutyl Salicylate is widely utilized in cosmetics and skincare formulations.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	5.1485 mL	25.7427 mL	51.4854 mL
5 mM	1.0297 mL	5.1485 mL	10.2971 mL
10 mM	0.5149 mL	2.5743 mL	5.1485 mL
50 mM	0.103 mL	0.5149 mL	1.0297 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Lapczynski A, et al. Fragrance material review on isobutyl salicylate. Food Chem Toxicol. 2007;45 Suppl 1:S424-7.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481