

Acetylcholinesterase Assay Kit (Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Acetylcholinesterase (AChE) belongs to serine hydrolases, which widely exists in various animal tissues and serum. AChE catalyzes the hydrolysis of acetylcholine (ACh) and plays an important role in the regulation of nerve conduction.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481