

## Acid Phosphatase Assay Kit (Spectrophotometry)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Acid Phosphatase (ACP) catalyzes the hydrolysis of phosphomonoesters into inorganic phosphate under acidic conditions. It is commonly found in the lysosomes of macrophages and is often used as an auxiliary diagnostic tool for prostate cancer.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481