

AH Assay Kit (Microanalysis)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Cytochrome P450 enzymes are a group of isozymes mainly present in the liver, which play an important role in the metabolism of exogenous substances, especially the metabolism of drugs and toxicants. Aniline-4-hydroxylase (AH) corresponds to the CYP2E1 subtype in the P450 enzyme system. CYP2E1 not only participates in drug metabolism but also catalyzes the activation process of various procarcinogens and protoxicants.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481