

Alcohol Dehydrogenase Assay Kit (Microanalysis)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Store at low temperature**
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|-------------|---|
| Description | Alcohol dehydrogenase (ADH) is a class of enzymes that catalyze the conversion of alcohols to aldehydes, accompanied by the reduction of NAD ⁺ to NADH. In humans, there are 9 types of ADH isozymes, with most of the ADH activity occurring in the liver. Members of the ADH family are the main enzymes involved in alcohol detoxification; genetic variations in ADH enzymes lead to differences in ADH activity and alcohol tolerance, and may regulate the susceptibility to alcoholism. |
|-------------|---|

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481