

Alkali Proteinase Assay Kit (Spectrophotometry)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Alkali Proteinase (AKP) refers to a class of enzymes that catalyze the hydrolysis of peptide bonds in proteins under alkaline conditions, belonging to serine proteases. In addition, this enzyme can also hydrolyze ester bonds and amide bonds, and has the functions of transesterification and transpeptidation. It is one of the main industrial enzymes, widely used in pharmaceutical, silk, food, leather-making and other industries.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481