

## Alkaline Phosphatase Assay Kit (Spectrophotometry)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Alkaline phosphatase (AKP/ALP) is a zinc-containing glycoprotein enzyme that can hydrolyze various natural and synthetic phospholipid monoester compounds in an alkaline environment. AKP/ALP is widely distributed in various organs of the human body, with the liver being the main one.
-------------	---

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481