

## Amylopectin Assay Kit (Microanalysis)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder -20°C for 3 years  
In solvent -80°C for 2 years

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Amylopectin is a highly branched polysaccharide. The ratio and content of amylose and amylopectin in starch have a direct impact on the processing, physicochemical properties, gelatinization temperature, etc. of starch products. Therefore, the research on starches with different ratios of amylose and amylopectin is of great significance.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481